

Vibrating Kegel Balls Operating Instructions



ATHENA'S
HOME NOVELTIES

SKU 3399

Designed by an all-women team, you'll appreciate every detail as you strengthen your vaginal walls and pelvic floor muscles! The two weighted spheres are covered in buttery soft silicone and weigh in at 2.25 ounces. Insert one or both and SQUEEZE! The flexible, 5.5" retrieval cord makes removable easy. Looking to add a WOW factor to your workout? They vibrate! 12 functions of vibration, escalation and pulsation will keep your mind on your muscle. Building strength has never been so pleasurable! The state-of-the-art memory chip knows what you like and resumes on the last function used! Your workouts can be wet and wild too! They are waterproof! Keep your vagina strong, healthy and happy with our *Rechargeable Kegel Balls!*

Features :

- Phthalate free
- 12 functions
- Memory chip resumes on last function used
- Total Length: 3.5in
- Larger Kegel Ball: 1.5in; Smaller Kegel Ball: 1.25in
- Silicone, ABS plastic detailing
- Rechargeable – USB cable included
- Not for use with silicone lubricants
- Waterproof



Cleaning and storage

Always clean product before first use, and clean after every use. For best results we recommend using Athena's (anti-microbial) Mighty Tidy Toy Cleaner (SKU MTT-1235 4 oz. or SKU MTT-2762 8 oz.)

Lubrication

Be sure to use ample lubrication with your bedroom toys. For great results try Athena's Perfect Harmony (SKU 2710 Original, SKU 2711 Tingling or SKU 2768 Warming). Not for use with silicone lubricants.

To ensure that you always enjoy a pleasurable experience with your adult novelty items, follow the care instructions listed below:

Do not store products near extreme cold or heat.

Store in a cool, dry location.

Store products separately as different materials may react with each other.

Take caution when placing lotion and massage oil bottles or adult novelties directly on furniture surfaces. Place them on a cloth to protect your toys and your furniture.

Warranty

All sales are final due to the intimate nature of our products. Motorized products are guaranteed for 30 days for a replacement. No refunds will be given. If you receive a broken product, or if your product breaks within 30 days of receipt, please contact our Customer Care Liaison for instructions on how to send back your item for a replacement. All defective product returns must be accompanied by a pre-authorization number. Our Customer Care Liaison can be reached by calling (401) 248-4721 or via email at E-mail: customerservice@athenashn.com



ATHENA'S
HOME NOVELTIES

www.athenashn.com

©September 2018 Athena's Home Novelties, Inc.
640 Winter Street Woonsocket RI 02895 1-877-ATHENAS

Charging

Charge upon opening. Connect the USB cable to the USB port of an AC Adapter, your computer, or other power source and the other end to the toys input port.

The red light will flash while charging and stays solid when fully charged. Approximately 45 minutes to fully charge. There is no danger of over-charging. Over 1 hour of playtime when on high speed.

Charge regularly every 6 months to preserve battery life.

Operating

To turn on or off, hold the power button down for 3 seconds.

Press the button to cycle through 12 intense functions.

Memory chip resumes on last function used.

Vibrating Kegel Ball will run for 30 minutes on highest speed, and 40 minutes on lowest speed.



**Charging Port
Located on Smaller Ball**

Power and Functions Button